medyglōbal

TEST
COVID-19
IgM/IgG Combo

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[Clinical evaluation]

Performance characteristic for the STANDARD Q COVID-19 IgM/IgG Combo Test for rapid detection of anti-SARS-CoV-2 antibodies was established in retrospective, single institute, randomized, single-blinded study conducted at a trial site in KOREA during the 2020 SARS-CoV-2 pandemic situation. A total of 336 retrospective specimens were tested using the STANDARD Q COVID-19 IgM/IgG Combo Test. These specimens consisted of serum from PCR positive or negative confirmed patients. The performance of the STANDARD Q COVID-19 IgM/IgG Combo Test were compared to a commercialized molecular assay. Although the STANDARD Q COVID-19 IgM/IgG Combo Test allows to test for IgM and IgG separately, due to the differing inter-patient time response to the virus, any individual with positive result for the IgM or the IgG test should be read as a positive for anti-SARS-CoV-2 antibodies. The combined test result (positive for IgM and/or IgG or negative for IgM and/or IgG) was used to calculate the total test sensitivity and specificity.

• Test sensitivity

The seroconversion time of IgM and IgG antibodies varies from person to person, but it was estimated to be around 7 days after onset of symptom^{4,5}. The STANDARD Q COVID-19 IgM/IgG Combo Test showed 96.50% of sensitivity using specimens from patients 7 days after symptom onset (combined IgM+IgG).

Table 1. Summary of the sensitivity of the STANDARD Q COVID-19 IgM/IgG Combo Test compared to PCR using specimens **7 days after symptom onset**

≥ 7 days after symptom onset		PCR		
		Positive	Negative	Total
STANDARD Q COVID-19 IgM/IgG Combo Test	Positive	138	0	138
	Negative	5	0	5
	Total	143	0	143
Sensitivity		96.50% (138/143, 95% CI, 92.03% - 98.86%)		

Table 2. Summary of the sensitivity of the STANDARD Q COVID-19 IgM/IgG Combo Test compared to PCR using specimens **within less than 7 days after symptom onset**

< 7 days after symptom onset		PCR		
		Positive	Negative	Total
STANDARD Q COVID-19 IgM/IgG Combo Test	Positive	26	0	26
	Negative	7	0	7
	Total	33	0	33
Sensitivity		78.79% (26/33, 95% CI, 61.09% - 91.02%)		

Table 3. Summary of the sensitivity of the STANDARD Q COVID-19 IgM/IgG Combo Test compared to PCR using specimens **in the period of 7 days to 14 days after symptom onset**

7-14 days after symptom onset		PCR		
		Positive	Negative	Total
STANDARD Q COVID-19 IgM/IgG Combo Test	Positive	52	0	52
	Negative	4	0	4
	Total	56	0	56
Sensitivity		92.86% (52/56, 95% CI, 82.71% - 98.02%)		

Table 4. Summary of the sensitivity of the STANDARD Q COVID-19 IgM/IgG Combo Test compared to PCR using specimens **14 days after symptom onset**

> 14 days after symptom onset		PCR		
		Positive	Negative	Total
STANDARD Q COVID-19 IgM/IgG Combo Test	Positive	86	0	86
	Negative	1	0	1
	Total	87	0	87
Sensitivity		98.85% (86/87, 95% CI, 93.76% - 99.97%)		

• Test specificity

The STANDARD Q COVID-19 IgM/IgG Combo Test showed **95% of specificity** (combined IgM+IgG).

Table 5. Summary of the specificity of the STANDARD Q COVID-19 IgM/IgG Combo Test compared to PCR

		PCR		
		Positive	Negative	Total
STANDARD Q COVID-19 IgM/IgG Combo Test	Positive	0	8	8
	Negative	0	152	152
	Total	0	160	160
Specificity		95.00% (152/160, 95% CI, 90.39% - 97.82%)		

ANALYTICAL PERFORMANCE

1. Limit of Detection: IgM-9.37 ug/ml, IgG-3.75 ug/ml

2. Cross-Reactivity: No cross-reactivity for the specimen of HIV Seroconversion panel, Japanese Encephalitis positive, Zika virus positive, Chikungunya positive, Dengue IgM positive, Salmonella typhi IgM positive, Rubella IgM, CMV IgG/IgM, Tick borne encephalitis IgM positive, West Nile Virus positive, Treponema palladium, HAV IgM positive, HAV IgG positive, HBV Ab positive, HCV Ab positive, Influenza vaccine positive, Leishmania positive, Brucella IgM positive, Chagas positive, Toxoplasma positive, Filariasis positive, Mycoplasma pneumonia IgM positive, Mycoplasma pneumonia IgG positive, Mycoplasma pneumonia IgM positive, Influenza A IgM positive, Influenza B IgM positive, Influenza A and B IgG+IgM positive and Tuberculosis positive.

3. Interference study: For Exogenous substances, there was no interference with medicines such as Zanamivir (for Influenza), Oseltamivir (for Influenza), Artemether-lumefantrine (for Malaria), Doxycycline hyclate (for Malaria), Quinine (for Malaria), Lamivudine (Retroviral medication), Ribavirin (for HCV), Daclatasvir (for HCV). Also there was no interference with Anti-inflammatory medication (Acetaminophen, Acetylsalicylic acid, Ibuprofen), Antibiotic (Erythromycin, Ciprofloxacin), Dietary substances (Caffeine, Ethanol, Biotin), Triglycerides, Cholesterol, Bilirubin (Unconjugate), Hemoglobin and Anticoagulants (EDTA, Heparin, Sodium citrate). For Endogenous substances, there was no interference with human anti-mouse antibody, whole blood of pregnant women and whole blood having elevated levels of hemoglobin, elevated levels of C-reactive protein and Elevated IgG or IgM samples.

4. Matrix Equivalency: The matrix and anticoagulant does not affect the detection of COVID-19 IgG and IgM in 25 contrived specimen between serum, plasma(Heparin, EDTA, Sodium citrate), venous whole blood(Heparin, EDTA, Sodium citrate), and capillary whole blood (collected in EDTA-treated tubes).

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1. From the Radiology Department of Sichuan Academy of Medical Science, Sichuan Provincial People's Hospital, Chengdu, China (W.K.); University of Michigan, Ann Arbor, MI, USA (P.P.A.).
2. Tahir ul Qamar, M.; Alqahtani, S.M.; Alamri, M.A.; Chen, L. Structural Basis of SARS-CoV-2 3CLpro and Anti-COVID-19 Drug Discovery from Medicinal Plants . Preprints 2020, 2020020193 (doi: 10.20944/preprints202002.0193.v1).

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Oficina central : C-4th&5th, 16, Deogyong-daero 1556beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16690, REPUBLICA DE COREA

Lugar de fabricación : 74, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, REPUBLICA DE COREA

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MT Promedt Consulting GmbH